

MEDLINE**UNITE**® Pre-Hydrated Bioimplant Technology Intelligently designed to maximize efficiency in flatfoot reconstruction



## Save Valuable O.R. Time.

In surgery, minutes matter. Time spent rehydrating and cutting freeze-dried bone is time wasted. Medline UNITE pre-hydrated bioimplant technology\* offers fully hydrated and pre-shaped wedges designed to enhance surgical speed, strength and performance.

## Pre-Hydrated for Speed and Strength

The patent-pending aqueous calcium chloride solution allows wedges to be stored fully hydrated for immediate use.

- » Eliminates idle time
- » Preserves structural integrity of the graft
- » Reduces the likelihood of intra- and post-operative graft crumbling

## Pre-Shaped for Uncompromised Performance

Made of dense cancellous bone, each wedge is pre-shaped to eliminate the time and waste of cutting a bone block.

- » Withstands the physical demands placed on structural grafts
- » Fully incorporates and resorbs—but removes easily if needed
- » Enables ancillary zero-profile screw fixation for rotational control

## Ordering Info

Item No.	Description	Size	Pkg
MWCT0005	Cotton Wedge	5mm	ea
MWCT0006		6mm	ea
MWCT0007		7mm	ea
MWEV0006	Evans Wedge	6mm	ea
MWEV0008		8mm	ea
MWEV0010		10mm	ea
MWET0012		12mm	ea
MWUT0012	Utility Wedge, (Tricortical)	12mm	ea

For more information or to schedule a case, contact your MEDLINEUNITE Representative or visit medlineunite.com

For use with:

Articulating Pin Distractor Allows you to adjust the correction, while providing unobstructed access to the osteotomy.

Wedge Trials

Allows you to view the correction visually and fluoroscopically before selecting the appropriate size wedge.

Evans and Cotton Plates

Available for ancillary fixation.



FOLLOW US fin & BLOG

MEDLINE

Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States 1-800-MEDLINE (633-5463) medline.com | info@medline.com Medline Canada 1-800-396-6996 medline.ca | canada@medline.com Medline México 01-800-831-0898 medlinemexico.com | mexico@medline.com