



TRAY LAYOUT GUIDE.

Ankle Fracture
Plating System

1 Lateral Fibula and Malleolar Hook Plates Caddy

Lateral Fibula Plates and Malleolar Hook Plates with 4.0mm Hook Plate Cannulated Screws

2 Posterior Fibula and Posterolateral Tibia Plates Caddy

Posterior Fibula and Posterolateral Tibia Plates

3 Straight Tibia-Fibula Plates Caddy

Straight Tibia-Fibula Plates and Flex Fibula Plates

4 Plate Screws Caddy

2.7mm Locking and Non-Locking Screws, 3.5mm Locking and Non-Locking Screws, 4.0mm Cancellous Screws with Temporary Fixation Pins

5 Instrument Middle Level


Bone Fragment Pick, Reduction Forceps, Ratchet AO/QC Handle, Pickups, Polyaxial Drill Guide, Compression Slot Drill Guides, Locking Drill Guides, Solid Core Drill Bits, Countersink, AO Style Depth Gauge, DrillGauge System, Straight and Retaining Drivers

6 Instrument Bottom Level


Bending Pliers, Bending Irons, Periosteal Elevator, Syndesmotomic Clamp, Malleolar Hook Plate Guide, Malleolar Hook Plate Impactor, Fibula Distractor, Distractor Pins

7 Cannulated Screws Caddy

4.0mm Cannulated Headed Screws with instruments (Guidewires, Cannulated Depth Gauge, Drill Bit, Driver, Washers)




Section 1



Lateral Fibula Plates

| Item No. | Description | Qty. |
|----------|--------------------|------|
| MPPA100L | Extra Small, Left | 2 |
| MPPA100R | Extra Small, Right | 2 |
| MPPA101L | Small, Left | 2 |
| MPPA101R | Small, Right | 2 |
| MPPA102L | Medium, Left | 2 |
| MPPA102R | Medium, Right | 2 |
| MPPA103L | Large, Left | 2 |
| MPPA103R | Large, Right | 2 |



Malleolar Hook Plates


| Item No. | Description | Qty. |
|----------|-----------------|------|
| MPPA301L | Small, Left | 2 |
| MPPA301R | Small, Right | 2 |
| MPPA302L | Standard, Left | 2 |
| MPPA302R | Standard, Right | 2 |



4.0mm, Cannulated Screws


| Item No. | Description | Qty. |
|----------|-------------|------|
| MPSH4040 | 4.0 X 40 mm | 2 |
| MPSH4045 | 4.0 X 45 mm | 2 |
| MPSH4050 | 4.0 X 50 mm | 2 |

Section 2



Posterior Fibula Plates


| Item No. | Description | Qty. |
|----------|--------------------|------|
| MPPA200L | Extra Small, Left | 1 |
| MPPA200R | Extra Small, Right | 1 |
| MPPA201L | Small, Left | 1 |
| MPPA201R | Small, Right | 1 |
| MPPA202L | Medium, Left | 2 |
| MPPA202R | Medium, Right | 2 |
| MPPA203L | Large, Left | 2 |
| MPPA203R | Large, Right | 2 |



Posterolateral Tibia Plates


| Item No. | Description | Qty. |
|----------|-----------------|------|
| MPPA401L | Small, Left | 2 |
| MPPA401R | Small, Right | 2 |
| MPPA402L | Standard, Left | 2 |
| MPPA402R | Standard, Right | 2 |

Section 3



Straight Tibia-Fibula Plates

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPPA503U | 3 Hole | 2 |
| MPPA504U | 4 Hole | 1 |
| MPPA505U | 5 Hole | 1 |
| MPPA506U | 6 Hole | 1 |
| MPPA507U | 7 Hole | 1 |
| MPPA508U | 8 Hole | 1 |
| MPPA509U | 9 Hole | 2 |
| MPPA510U | 10 Hole | 2 |
| MPPA511U | 11 Hole | 2 |
| MPPA512U | 12 Hole | 2 |



Flex Fibula Plates

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPPA004U | 4 Hole | 1 |
| MPPA005U | 5 Hole | 1 |
| MPPA006U | 6 Hole | 1 |
| MPPA007U | 7 Hole | 1 |
| MPPA008U | 8 Hole | 1 |

Section 4




2.7 mm, Polyaxial Locking Screws

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPSL2710 | 2.7 X 10 mm | 4 |
| MPSL2712 | 2.7 X 12 mm | 4 |
| MPSL2714 | 2.7 X 14 mm | 4 |
| MPSL2716 | 2.7 X 16 mm | 4 |
| MPSL2718 | 2.7 X 18 mm | 4 |
| MPSL2720 | 2.7 X 20 mm | 4 |




2.7 mm, Non-Locking Screws

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPSN2710 | 2.7 X 10 mm | 2 |
| MPSN2712 | 2.7 X 12 mm | 2 |
| MPSN2714 | 2.7 X 14 mm | 2 |
| MPSN2716 | 2.7 X 16 mm | 2 |
| MPSN2718 | 2.7 X 18 mm | 2 |
| MPSN2720 | 2.7 X 20 mm | 2 |



3.5 mm, Polyaxial Locking Screws

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPSL3510 | 3.5 X 10 mm | 4 |
| MPSL3512 | 3.5 X 12 mm | 4 |
| MPSL3514 | 3.5 X 14 mm | 4 |
| MPSL3516 | 3.5 X 16 mm | 4 |
| MPSL3518 | 3.5 X 18 mm | 4 |
| MPSL3520 | 3.5 X 20 mm | 4 |
| MPSL3522 | 3.5 X 22 mm | 4 |
| MPSL3524 | 3.5 X 24 mm | 4 |
| MPSL3526 | 3.5 X 26 mm | 4 |
| MPSL3528 | 3.5 X 28 mm | 4 |
| MPSL3530 | 3.5 X 30 mm | 4 |
| MPSL3532 | 3.5 X 32 mm | 4 |
| MPSL3534 | 3.5 X 34 mm | 4 |
| MPSL3536 | 3.5 X 36 mm | 4 |
| MPSL3538 | 3.5 X 38 mm | 4 |
| MPSL3540 | 3.5 X 40 mm | 4 |
| MPSL3542 | 3.5 X 42 mm | 2 |
| MPSL3544 | 3.5 X 44 mm | 2 |
| MPSL3546 | 3.5 X 46 mm | 2 |
| MPSL3548 | 3.5 X 48 mm | 2 |
| MPSL3550 | 3.5 X 50 mm | 2 |
| MPSL3555 | 3.5 X 55 mm | 2 |
| MPSL3560 | 3.5 X 60 mm | 2 |



3.5 mm, Non-Locking Screws

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPSN3510 | 3.5 X 10 mm | 6 |
| MPSN3512 | 3.5 X 12 mm | 6 |
| MPSN3514 | 3.5 X 14 mm | 6 |
| MPSN3516 | 3.5 X 16 mm | 6 |
| MPSN3518 | 3.5 X 18 mm | 6 |
| MPSN3520 | 3.5 X 20 mm | 6 |
| MPSN3522 | 3.5 X 22 mm | 4 |
| MPSN3524 | 3.5 X 24 mm | 4 |
| MPSN3526 | 3.5 X 26 mm | 4 |
| MPSN3528 | 3.5 X 28 mm | 4 |
| MPSN3530 | 3.5 X 30 mm | 4 |
| MPSN3532 | 3.5 X 32 mm | 4 |
| MPSN3534 | 3.5 X 34 mm | 4 |
| MPSN3536 | 3.5 X 36 mm | 4 |
| MPSN3538 | 3.5 X 38 mm | 4 |
| MPSN3540 | 3.5 X 40 mm | 4 |
| MPSN3542 | 3.5 X 42 mm | 2 |
| MPSN3544 | 3.5 X 44 mm | 2 |
| MPSN3546 | 3.5 X 46 mm | 2 |
| MPSN3548 | 3.5 X 48 mm | 2 |
| MPSN3550 | 3.5 X 50 mm | 2 |
| MPSN3552 | 3.5 X 52 mm | 2 |
| MPSN3554 | 3.5 X 54 mm | 2 |
| MPSN3556 | 3.5 X 56 mm | 2 |
| MPSN3558 | 3.5 X 58 mm | 2 |
| MPSN3560 | 3.5 X 60 mm | 2 |

Section 4



4.0 mm, Cancellous Screws

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPSC4010 | 4.0 X 10 mm | 4 |
| MPSC4012 | 4.0 X 12 mm | 4 |
| MPSC4014 | 4.0 X 14 mm | 4 |
| MPSC4016 | 4.0 X 16 mm | 4 |
| MPSC4018 | 4.0 X 18 mm | 4 |
| MPSC4020 | 4.0 X 20 mm | 4 |



Temporary Fixation Pins

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPPF1115 | 1.1 X 15 mm | 6 |

Section 5



Bone Fragment Pick

| Item No. | Description | Qty. |
|----------|--------------------|------|
| MPN50005 | Bone Fragment Pick | 1 |




Ratchet AO/QC Handle

| Item No. | Description | Qty. |
|----------|-------------|------|
| MSN90001 | Cannulated | 1 |



Polyaxial Drill Guide

| Item No. | Description | Qty. |
|----------|----------------|------|
| MPN40007 | ● 2.0 ● 2.8 mm | 1 |



Locking Tower Drill Guides

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPN40001 | ● 2.0 mm | 2 |
| MPN40002 | ● 2.8 mm | 2 |



Reduction Forceps

| Item No. | Description | Qty. |
|----------|---------------------|------|
| MPN50001 | Pointed, Osteopenic | 1 |
| MPN50002 | Pointed, Wide | 1 |
| MPN50003 | Serrated | 2 |
| MPN50015 | Pointed | 1 |



Pickups

| Item No. | Description | Qty. |
|----------|-------------|------|
| MSN90003 | Pickups | 1 |

Section 5




Compression Slot Drill Guides

| Item No. | Description | Qty. |
|----------|-----------------|------|
| MPN50011 | ● 2.0 mm/2.8 mm | 1 |
| MPN50012 | ● 2.8 mm/3.5 mm | 1 |



Countersink


| Item No. | Description | Qty. |
|----------|------------------|------|
| MSN20003 | ● 3.0/3.5 Heated | 1 |



DrillGauge 3-in-1 System

| Item No. | Description | Qty. |
|----------|--------------------|------|
| MPN50010 | ● Assembly, 2.8 mm | 1 |
| MPN10128 | ● Drill Bit, 2.8mm | 2 |

*A 2.0mm DrillGauge 3-in-1 System is also available upon request.




Solid Core Drill Bits

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPN10020 | ● 2.0 mm | 2 |
| MPN10028 | ● 2.8 mm | 2 |
| MPN10035 | ● 3.5 mm | 2 |



AO Style Depth Gauge


| Item No. | Description | Qty. |
|----------|-------------|------|
| MPN40006 | 60 mm | 1 |



Drivers

| Item No. | Description | Qty. |
|----------|------------------|------|
| MPN30001 | ● Straight, T15 | 2 |
| MPN30002 | ● Retaining, T15 | 2 |

Section 6



Bending Pliers

| Item No. | Description | Qty. |
|----------|----------------|------|
| MPN50008 | Bending Pliers | 2 |



Bending Irons

| Item No. | Description | Qty. |
|----------|---------------|------|
| MPN50007 | Bending Irons | 2 |



Periosteal Elevator

| Item No. | Description | Qty. |
|----------|---------------------|------|
| MPN50006 | Periosteal Elevator | 1 |



Syndesmotomic Clamp

| Item No. | Description | Qty. |
|----------|-----------------|------|
| MPN50004 | Malleolar Clamp | 1 |

Section 6



Malleolar Hook Plate Guide

| Item No. | Description | Qty. |
|----------|------------------|------|
| MPN50014 | Hook Plate Guide | 1 |



Malleolar Hook Plate Impactor

| Item No. | Description | Qty. |
|----------|---------------------|------|
| MPN50009 | Hook Plate Impactor | 1 |



Distractor Pins

| Item No. | Description | Qty. |
|----------|--------------|------|
| MPX27100 | 2.7 X 100 mm | 4 |



Fibula Distractor

| Item No. | Description | Qty. |
|----------|-------------|------|
| MPX00002 | Distractor | 1 |

Section 7



4.0 mm, Cannulated Headed Screws

| Item No. | Description | Qty. |
|----------|-------------|------|
| MSD04030 | 4.0 X 30 mm | 2 |
| MSD04032 | 4.0 X 32 mm | 2 |
| MSD04034 | 4.0 X 34 mm | 2 |
| MSD04036 | 4.0 X 36 mm | 2 |
| MSD04038 | 4.0 X 38 mm | 2 |
| MSD04040 | 4.0 X 40 mm | 2 |
| MSD04042 | 4.0 X 42 mm | 2 |
| MSD04044 | 4.0 X 44 mm | 2 |
| MSD04046 | 4.0 X 46 mm | 2 |
| MSD04048 | 4.0 X 48 mm | 2 |
| MSD04050 | 4.0 X 50 mm | 2 |
| MSD04055 | 4.0 X 55 mm | 2 |
| MSD04060 | 4.0 X 60 mm | 2 |



Cannulated Drill Bit

| Item No. | Description | Qty. |
|----------|----------------------|------|
| MSN10011 | ●● Hook Screw, 3.0mm | 2 |



Washers

| Item No. | Description | Qty. |
|----------|------------------------------|------|
| MSW04045 | For 4.0/4.5 mm Headed Screws | 4 |

Non-Threaded Guidewires

| Item No. | Description | Qty. |
|----------|--------------|------|
| MSG14150 | 1.4 X 150 mm | 6 |



Cannulated Depth Gauge

| Item No. | Description | Qty. |
|----------|-------------|------|
| MSN40001 | 150 mm Wire | 1 |



Cannulated Driver

| Item No. | Description | Qty. |
|----------|--------------------|------|
| MSN30003 | ●● Cannulated, T20 | 2 |



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093
 Medline United States
 1-800-MEDLINE (633-5463)
 medlineunite.com | unite@medline.com

Medline Canada
 1-800-396-6996 | medline.ca | canada@medline.com
Medline México
 01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    