

FIXATED ON PRECISION.

The difference is in the details.

We know that small details make a big difference in driving precision and efficiency. That's why our snap-off screws are so meticulously crafted. We thoughtfully fine-tune every element to help improve outcomes and elevate performance.

MEDLINE**UNITE**® Snap-Off Screws

- » Self-cutting feature under head for positioning flush with proximal cortex, even at an oblique angle
- » Combined utility of countersinking driver and self-cutting screws reduces likelihood of soft tissue irritation and helps improve outcomes
- » Ideal for fixation of Weil metatarsal osteotomies, bunionectomies and small bone osteotomies
- » Available in 2.0mm and 2.7mm; color-coded to for easy identification

Item No.	Description
SNAP2011	Snap-Off Screw, 2.0 x 11 mm
SNAP2012	Snap-Off Screw, 2.0 x 12 mm
SNAP2013	Snap-Off Screw, 2.0 x 13 mm
SNAP2014	Snap-Off Screw, 2.0 x 14 mm
SNAP2713	Snap-Off Screw, 2.7 x 13 mm
SNAP2715	Snap-Off Screw, 2.7 x 15 mm
SNAP2717	Snap-Off Screw, 2.7 x 17 mm
SNAP2719	Snap-Off Screw, 2.7 x 19 mm
SNAP2721	Snap-Off Screw, 2.7 x 21 mm
SNAPWIRE	Guidewire, 1.1 x 50 mm
SNAPDRSS	Driver 3-in-1 Instrument, AO/QC
SNAPHAND	Driver 3-in-1 Instrument, Handle

For more information or to schedule a case, contact your

MEDLINEUNITE Representative or visit medlineunite.com



Compact Caddy

Sturdy, yet low profile with a small footprint. Sterilize and use independently, or house within MEDLINE**UNITE** Cannulated Screw Trays if multiple osteotomies are being performed.



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States 1-800-MEDLINE (633-5463) medline.com | info@medline.com Medline Canada 1-800-396-6996 medline.ca | canada@medline.com Medline México 01-800-831-0898 medlinemexico.com | mexico@medline.com

FOLLOW US fin | BLOG