



UNITE[®]
FOOT & ANKLE

Visit us at
AOFAS
Hyatt Regency Chicago
Booth #229
September 12-15

INTELLIGENTLY DESIGNED

DRILLED INTO ACCURACY.

Lateral Fibula Plate with
Fibula Lengthening Distractor
and DrillGauge




DrillGauge 3-in-1 spring-loaded system acts as a drill guide, drill bit and depth gauge for faster, one-handed technique.



Distractor locks into plate for easier, more controlled anatomic reduction—modified “push-pull” technique.


PART OF THE UNITE ANKLE FRACTURE PLATING SYSTEM.



LATERAL FIBULA



FLEX FIBULA



STRAIGHT TIBIA-FIBULA



POSTERIOR FIBULA,
POSTEROLATERAL TIBIA
AND MALLEOLAR HOOK

For more information or to schedule a case, contact your Medline UNITE representative, or visit medlineunite.com.