



INSTRUMENTED INNOVATION




Distal Tibia Plate Inserter features built-in locking drill guides allowing for MIPO technique when Pilon fractures extend more proximally.



Fibula Lengthening Distractor locks into place for easier, more controlled anatomic reduction – modified “push-pull” technique.



Sinus Tarsi Extension Plate Inserter allows for easy plate insertion and positioning through a sinus tarsi incision for percutaneous screw placement in the posterior tuberosity.



UNITE[®]
FOOT & ANKLE



Scan this code to learn more about
the Ankle Fracture Plating System.