



# TRAY LAYOUT GUIDE.

Ankle Fusion  
Plating System

## 1 Ankle Fusion Plate Caddy

Anterior Plates, Posterior and Lateral Plates, Targeting Guide, AO/QC Non-Ratcheting Handle, Curettes and Arthrotomes

## 2 Plate Screw Caddy

Locking and Non-Locking Screws (3.5mm, 4.0mm, 4.5mm, 5.5mm), Temporary Fixation Pins, Locking Drill Guides, and Retaining Drivers

## 3 Ankle Distraction


Two ankle joint distractors and 2.5mm x 100mm distraction pins

## 4 Instruments Bottom Level

Hohmann Retractors, Tissue Protector Drill Guide/Poly Axial Cones, Compression Slot Drill Guide, Drills, Depth Gauge, AO/QC Ratcheting Handle, Pick-ups, Cross-Plate Drill Guide, Fenestration Drill Pin

## 5 Extra Long Implants Caddy

Primary Pilon Plates, 6.5mm Headless Cannulated Screws, 7.0mm Headed Cannulated Screws



# Section 1



## Standard Anterior Plates

Item No.	Description	Qty.
MPAF201L	Small, Left	1
MPAF201R	Small, Right	1
MPAF202L	Standard, Left	1
MPAF202R	Standard, Right	1



## Petite Inline Plates

Item No.	Description	Qty.
MPAF001U	Small	1
MPAF002U	Standard	1



## Petite T-Style Plates

Item No.	Description	Qty.
MPAF003U	Small	1
MPAF004U	Standard	1



## Anterior Short Talar Neck Plates

Item No.	Description	Qty.
MPAF101L	Small, Left	1
MPAF101R	Small, Right	1
MPAF102L	Standard, Left	1
MPAF102R	Standard, Right	1

# Section 1




## Posterior TT/TTC Plates

Item No.	Description	Qty.
MPAF301L	TT, Left	1
MPAF301R	TT, Right	1
MPAF401L	TTC, Left	1
MPAF401R	TTC, Right	1



## Lateral Plates

Item No.	Description	Qty.
MPAF501L	Standard, Left	1
MPAF501R	Standard, Right	1



## Targeting Guide & Guidepin Sleeves

Item No.	Description	Qty.
MPN52015	Targeting Guide Arm, Left	1
MPN52017	Targeting Guide Arm, Right	1
MPN52016	Targeting Guide Sleeve, 2.0mm	1
MPN52018	Targeting Guide Sleeve, 2.5mm	1

# Section 1



## AO/QC Non-Ratcheting Handle

Item No.	Description	Qty.
MSN90006	AO/QC Non-Ratcheting Cannulated Handle	1



## Curettes

Item No.	Description	Qty.
MPN52010	Curved Cup Curette	1
MPN52011	Straight Cup Curette	1



## Arthrotomes

Item No.	Description	Qty.
MPN52012	6.5mm Arthrotome	1
MPN52013	13mm Arthrotome	1

# Section 2



## Temporary Fixation Pins

Item No.	Description	Qty.
MPPF2020	2.0 X 20mm, Smooth	4
MPPF2020T	2.0 X 20mm, Threaded	4



## Locking Drill Guides


Item No.	Description	Qty.
MPN40002	2.8mm	2
MPN40008	3.0mm	2
MPN52004	3.1mm	2
MPN52005	3.8mm	2



## Drivers


Item No.	Description	Qty.
MPN30002	Retaining, T15	2
MPN30005	Retaining, T20	2

# Section 2



## 3.5 mm, Polyaxial Locking Screws

Item No.	Description	Qty.
MPSL3516	3.5 X 16 mm	2
MPSL3518	3.5 X 18 mm	2
MPSL3520	3.5 X 20 mm	2
MPSL3522	3.5 X 22 mm	2
MPSL3524	3.5 X 24 mm	2
MPSL3526	3.5 X 26 mm	2
MPSL3528	3.5 X 28 mm	2
MPSL3530	3.5 X 30 mm	2
MPSL3532	3.5 X 32 mm	2
MPSL3534	3.5 X 34 mm	2
MPSL3536	3.5 X 36 mm	2
MPSL3538	3.5 X 38 mm	2
MPSL3540	3.5 X 40 mm	2



## 3.5 mm, Non-Locking Screws

Item No.	Description	Qty.
MPSN3516	3.5 X 16 mm	2
MPSN3518	3.5 X 18 mm	2
MPSN3520	3.5 X 20 mm	2
MPSN3522	3.5 X 22 mm	2
MPSN3524	3.5 X 24 mm	3
MPSN3526	3.5 X 26 mm	3
MPSN3528	3.5 X 28 mm	3
MPSN3530	3.5 X 30 mm	3
MPSN3532	3.5 X 32 mm	3
MPSN3534	3.5 X 34 mm	3
MPSN3536	3.5 X 36 mm	3
MPSN3538	3.5 X 38 mm	2
MPSN3540	3.5 X 40 mm	2



## 4.0mm Polyaxial Locking Screws

Item No.	Description	Qty.
MPSL4016	4.0 x 16mm	2
MPSL4018	4.0 x 18mm	2
MPSL4020	4.0 x 20mm	2
MPSL4022	4.0 x 22mm	2
MPSL4024	4.0 x 24mm	2
MPSL4026	4.0 x 26mm	2
MPSL4028	4.0 x 28mm	2
MPSL4030	4.0 x 30mm	2
MPSL4032	4.0 x 32mm	2
MPSL4034	4.0 x 34mm	2
MPSL4036	4.0 x 36mm	2
MPSL4038	4.0 x 38mm	2
MPSL4040	4.0 x 40mm	2



## 4.0mm Non-Locking Screws

Item No.	Description	Qty.
MPSN4016	4.0 x 16mm	2
MPSN4018	4.0 x 18mm	2
MPSN4020	4.0 x 20mm	2
MPSN4022	4.0 x 22mm	2
MPSN4024	4.0 x 24mm	2
MPSN4026	4.0 x 26mm	2
MPSN4028	4.0 x 28mm	2
MPSN4030	4.0 x 30mm	2
MPSN4032	4.0 x 32mm	2
MPSN4034	4.0 x 34mm	2
MPSN4036	4.0 x 36mm	2
MPSN4038	4.0 x 38mm	2
MPSN4040	4.0 x 40mm	2

# Section 2



## 4.5mm Polyaxial Locking Screws

Item No.	Description	Qty.
MPSL4516	4.5 x 16mm	3
MPSL4518	4.5 x 18mm	3
MPSL4520	4.5 x 20mm	3
MPSL4522	4.5 x 22mm	3
MPSL4524	4.5 x 24mm	3
MPSL4526	4.5 x 26mm	3
MPSL4528	4.5 x 28mm	3
MPSL4530	4.5 x 30mm	3
MPSL4532	4.5 x 32mm	3
MPSL4534	4.5 x 34mm	3
MPSL4536	4.5 x 36mm	3
MPSL4538	4.5 x 38mm	3
MPSL4540	4.5 x 40mm	3
MPSL4542	4.5 x 42mm	2
MPSL4544	4.5 x 44mm	2
MPSL4546	4.5 x 46mm	2
MPSL4548	4.5 x 48mm	2
MPSL4550	4.5 x 50mm	2
MPSL4552	4.5 x 52mm	2
MPSL4554	4.5 x 54mm	2
MPSL4556	4.5 x 56mm	2
MPSL4558	4.5 x 58mm	2
MPSL4560	4.5 x 60mm	2



## 4.5mm Non-Locking Screws

Item No.	Description	Qty.
MPSN4516	4.5 x 16mm	2
MPSN4518	4.5 x 18mm	2
MPSN4520	4.5 x 20mm	2
MPSN4522	4.5 x 22mm	2
MPSN4524	4.5 x 24mm	2
MPSN4526	4.5 x 26mm	2
MPSN4528	4.5 x 28mm	2
MPSN4530	4.5 x 30mm	2
MPSN4532	4.5 x 32mm	2
MPSN4534	4.5 x 34mm	2
MPSN4536	4.5 x 36mm	2
MPSN4538	4.5 x 38mm	2
MPSN4540	4.5 x 40mm	2
MPSN4542	4.5 x 42mm	2
MPSN4544	4.5 x 44mm	2
MPSN4546	4.5 x 46mm	2
MPSN4548	4.5 x 48mm	2
MPSN4550	4.5 x 50mm	2
MPSN4552	4.5 x 52mm	2
MPSN4554	4.5 x 54mm	2
MPSN4556	4.5 x 56mm	2
MPSN4558	4.5 x 58mm	2
MPSN4560	4.5 x 60mm	2

# Section 2



## 5.5mm Polyaxial Locking Screws

Item No.	Description	Qty.
MPSL5516	5.5 x 16mm	3
MPSL5518	5.5 x 18mm	3
MPSL5520	5.5 x 20mm	3
MPSL5522	5.5 x 22mm	3
MPSL5524	5.5 x 24mm	3
MPSL5526	5.5 x 26mm	3
MPSL5528	5.5 x 28mm	3
MPSL5530	5.5 x 30mm	3
MPSL5532	5.5 x 32mm	3
MPSL5534	5.5 x 34mm	3
MPSL5536	5.5 x 36mm	3
MPSL5538	5.5 x 38mm	3
MPSL5540	5.5 x 40mm	3
MPSL5542	5.5 x 42mm	2
MPSL5544	5.5 x 44mm	2
MPSL5546	5.5 x 46mm	2
MPSL5548	5.5 x 48mm	2
MPSL5550	5.5 x 50mm	2
MPSL5552	5.5 x 52mm	2
MPSL5554	5.5 x 54mm	2
MPSL5556	5.5 x 56mm	2
MPSL5558	5.5 x 58mm	2
MPSL5560	5.5 x 60mm	2



## 5.5mm Non-Locking Screws

Item No.	Description	Qty.
MPSN5516	5.5 x 16mm	2
MPSN5518	5.5 x 18mm	2
MPSN5520	5.5 x 20mm	2
MPSN5522	5.5 x 22mm	2
MPSN5524	5.5 x 24mm	2
MPSN5526	5.5 x 26mm	2
MPSN5528	5.5 x 28mm	2
MPSN5530	5.5 x 30mm	2
MPSN5532	5.5 x 32mm	2
MPSN5534	5.5 x 34mm	2
MPSN5536	5.5 x 36mm	2
MPSN5538	5.5 x 38mm	2
MPSN5540	5.5 x 40mm	2
MPSN5542	5.5 x 42mm	2
MPSN5544	5.5 x 44mm	2
MPSN5546	5.5 x 46mm	2
MPSN5548	5.5 x 48mm	2
MPSN5550	5.5 x 50mm	2
MPSN5552	5.5 x 52mm	2
MPSN5554	5.5 x 54mm	2
MPSN5556	5.5 x 56mm	2
MPSN5558	5.5 x 58mm	2
MPSN5560	5.5 x 60mm	2



# Section 3



## Distractor

Item No.	Description	Qty.
MPX00003	Distractor	2



## Distractor Pin

Item No.	Description	Qty.
MPX25100	2.5 X 100mm	6

# Section 4



## Hohmann

Item No.	Description	Qty.
MPN52009	Hohmann Retractor	2



## Tissue Protector Drill Guide / Polyaxial Cone

Item No.	Description	Qty.
MPN52014	2.8mm	1
MPN52008	3.0mm	1
MPN52002	3.1mm	1
MPN52003	3.8mm	1



## Compression Slot Drill Guide

Item No.	Description	Qty.
MPN52007	3.1/3.8mm	1



## Fenestration Drill Pin

Item No.	Description	Qty.
MPDP0025	2.5mm	2



## Drill Bit

Item No.	Description	Qty.
MPN10028	2.8mm	2
MPN10030	3.0mm	2
MPN10031	3.1mm	2
MPN10038	3.8mm	2
MPN10045	4.5 mm	2



## Depth Gauge

Item No.	Description	Qty.
MPN52001	Depth Gauge	1

# Section 4



Ratcheting AO/QC Handle

Item No.	Description	Qty.
MSN90001	Cannulated	1



Pickups

Item No.	Description	Qty.
MSN90003	Pickups	1



Cross-Plate Drill Guide

Item No.	Description	Qty.
MPN52019	3.1mm/4.5mm	1

# Section 5



## Primary Pilon Fusion Plates

Item No.	Description	Qty.
MPAF601L	Left	1
MPAF601R	Right	1



## 6.5mm Cannulated Headless Screws

Item No.	Description	Qty.
MSL06505	6.5mm X 105mm	2
MSL06510	6.5mm X 110mm	2
MSL06515	6.5mm X 115mm	2
MSL06520	6.5mm X 120mm	2
MSL06525	6.5mm X 125mm	2



## 7.0mm Cannulated Headed Screws

Item No.	Description	Qty.
MSD07005	7.0mm X 105mm	2
MSD07010	7.0mm X 110mm	2
MSD07015	7.0mm X 115mm	2
MSD07020	7.0mm X 120mm	2
MSD07025	7.0mm X 125mm	2



**Medline Industries, LP.**  
 Three Lakes Drive, Northfield, IL 60093  
 Medline United States | 1-800-MEDLINE (633-5463)  
 medline.com | info@medline.com

**Medline Canada**  
 1-800-268-2848 | medline.ca | canada@medline.com  
**Medline México**  
 01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    