



TRAY LAYOUT GUIDE

Mini Foot Plating System

1 Mini Foot Plates Caddy

Mini Foot Plates (L Plates, T Plates, Straight Plates, Z Plates, Spiral Plates, Talar Neck, Cuboid, Navicular)

2 5th Met Hook Plate Caddy

5th Met Hook Plate System Plates with instruments (Nitinol Drill Pins, Guide, Inserter, Guidewires, Depth Gauge, Drill Bit, Driver)

3 Mini Plate Screws Caddy

2.0mm, 2.4mm, 2.7mm Locking and Non-Locking Mini Plate Screws, 4.0mm Cannulated Screw for Hook Plate, Temporary Fixation Pins

4 Instrument Rack Top Level


Locking Drill Guides, Solid Core Drill Bits, AO Style Depth Gauge, Ratchet AO/QC Handle, Pickups, Retaining Drivers

5 Instrument Rack Bottom Level (Left)


Bone Fragment Pick, Retractors, Reduction Forceps, Compression Slot Drill Guides, Polyaxial Drill Guide

6 Instrument Rack Bottom Level (Right)

Periosteal Elevator, 5th Met Hook Plate Impactor, Bending Pliers, Plate Cutter




Section 1




Metatarsal L Plate

Item No.	Description	Qty.
MPPM201U	4 Hole	2
MPPM202U	5 Hole	2
MPPM203U	7 Hole	2
MPPM204U	9 Hole	2



Metatarsal Z Plate

Item No.	Description	Qty.
MPPM400U	Small	2
MPPM401U	Medium	2
MPPM402U	Large	2



Metatarsal T Plate

Item No.	Description	Qty.
MPPM301U	4 Hole	3
MPPM302U	5 Hole	3
MPPM303U	7 Hole	3
MPPM304U	9 Hole	3




Talar Neck Plate

Item No.	Description	Qty.
MPPM800U	Universal	2



Cuboid Plate

Item No.	Description	Qty.
MPPM600U	Universal	2



Metatarsal Straight Plate

Item No.	Description	Qty.
MPPM104U	4 Hole	3
MPPM105U	5 Hole	3
MPPM107U	7 Hole	3
MPPM109U	9 Hole	3



Navicular Plate


Item No.	Description	Qty.
MPPM700L	Left	2
MPPM700R	Right	2

Section 2



5th Metatarsal Spiral Plate

Item No.	Description	Qty.
MPPM510L	Left	1
MPPM510R	Right	1



5th Metatarsal Hook Plate

Item No.	Description	Qty.
MPPM501L	Standard, Left	2
MPPM501R	Standard, Right	2
MPPM502L	Large, Left	2
MPPM502R	Large, Right	2



5th Met Hook Plate Guide

Item No.	Description	Qty.
MPN51014	Left	1
MPN51013	Right	1



5th Met Nitinol Drill Pins

Item No.	Description	Qty.
MPDP0017	1.7 mm	2

Non-Threaded Guidewires

Item No.	Description	Qty.
MSG14150	1.4 X 150 mm	6



Cannulated Depth Gauge

Item No.	Description	Qty.
MSN40001	150 mm Wire	1



Cannulated Drill Bit

Item No.	Description	Qty.
MSN10012	● Cannulated Screw, 2.6 mm	1



Cannulated Driver

Item No.	Description	Qty.
MSN30005	● Cannulated T15	1



5th Met Hook Plate Inserter

Item No.	Description	Qty.
MPN51009	Left	1
MPN51010	Right	1

Section 3



2.0 mm, Mini Polyaxial Locking Screws

Item No.	Description	Qty.
MMSL2010	2.0 x 10 mm	4
MMSL2012	2.0 x 12 mm	4
MMSL2014	2.0 x 14 mm	4
MMSL2016	2.0 x 16 mm	4
MMSL2018	2.0 x 18 mm	4
MMSL2020	2.0 x 20 mm	4



2.0 mm, Mini Non-Locking Screws

Item No.	Description	Qty.
MMSN2010	2.0 x 10 mm	4
MMSN2012	2.0 x 12 mm	4
MMSN2014	2.0 x 14 mm	4
MMSN2016	2.0 x 16 mm	4
MMSN2018	2.0 x 18 mm	4
MMSN2020	2.0 x 20 mm	4



2.4 mm, Mini Polyaxial Locking Screws

Item No.	Description	Qty.
MMSL2408	2.4 x 8 mm	6
MMSLL2410	2.4 x 10 mm	6
MMSL2412	2.4 x 12 mm	6
MMSL2414	2.4 x 14 mm	6
MMSL2416	2.4 x 16 mm	6
MMSL2418	2.4 x 18 mm	6
MMSL2420	2.4 x 20 mm	6
MMSL2422	2.4 x 22 mm	6
MMSL2424	2.4 x 24 mm	4
MMSL2426	2.4 x 26 mm	4
MMSL2428	2.4 x 28 mm	4
MMSL2430	2.4 x 30 mm	4



2.4 mm, Mini Non-Locking Screws

Item No.	Description	Qty.
MMSN2408	2.4 x 8 mm	4
MMSN2410	2.4 x 10 mm	4
MMSN2412	2.4 x 12 mm	4
MMSN2414	2.4 x 14 mm	4
MMSN2416	2.4 x 16 mm	4
MMSN2418	2.4 x 18 mm	4
MMSN2420	2.4 x 20 mm	4
MMSN2422	2.4 x 22 mm	2
MMSN2424	2.4 x 24 mm	2
MMSN2426	2.4 x 26 mm	2
MMSN2428	2.4 x 28 mm	2
MMSN2430	2.4 x 30 mm	2



2.7 mm, Mini Polyaxial Locking Screws

Item No.	Description	Qty.
MMSL2708	2.7 x 8 mm	4
MMSL2710	2.7 x 10 mm	4
MMSL2712	2.7 x 12 mm	4
MMSL2714	2.7 x 14 mm	4
MMSL2716	2.7 x 16 mm	4
MMSL2718	2.7 x 18 mm	4
MMSL2720	2.7 x 20 mm	4
MMSL2722	2.7 x 22 mm	4
MMSL2724	2.7 x 24 mm	4
MMSL2726	2.7 x 26 mm	4
MMSL2728	2.7 x 28 mm	4
MMSL2730	2.7 x 30 mm	4
MMSL2732	2.7 x 32 mm	4
MMSL2734	2.7 x 34 mm	4
MMSL2736	2.7 x 36 mm	4
MMSL2738	2.7 x 38 mm	2
MMSL2740	2.7 x 40 mm	2



2.7 mm, Mini Non-Locking Screws

Item No.	Description	Qty.
MMSN2708	2.7 x 8 mm	4
MMSN2710	2.7 x 10 mm	4
MMSN2712	2.7 x 12 mm	4
MMSN2714	2.7 x 14 mm	4
MMSN2716	2.7 x 16 mm	4
MMSN2718	2.7 x 18 mm	4
MMSN2720	2.7 x 20 mm	4
MMSN2722	2.7 x 22 mm	2
MMSN2724	2.7 x 24 mm	2
MMSN2726	2.7 x 26 mm	2
MMSN2728	2.7 x 28 mm	2
MMSN2730	2.7 x 30 mm	2
MMSN2732	2.7 x 32 mm	2
MMSN2734	2.7 x 34 mm	2
MMSN2736	2.7 x 36 mm	2
MMSN2786	2.7 x 86 mm	2
MMSN2740	2.7 x 40 mm	2



4.0 mm, Cannulated Hook Plate Screws

Item No.	Description	Qty.
MDS14020	4.0 x 20 mm	2
MDS14022	4.0 x 22 mm	2
MDS14024	4.0 x 24 mm	2
MDS14026	4.0 x 26 mm	2
MDS14028	4.0 x 28 mm	2
MDS14030	4.0 x 30 mm	2
MDS14032	4.0 x 32 mm	2
MDS14034	4.0 x 34 mm	2
MDS14036	4.0 x 36 mm	2



Temporary Fixation Pins

Item No.	Description	Qty.
MMPF1110	Smooth Mini, 1.1 x 10 mm	4
MMPF1110T	Threaded Mini, 1.1 x 10 mm	4

Section 4



Locking Tower Drill Guides

Item No.	Description	Qty.
MPN42101	1.6 mm	2
MPN42102	1.8 mm	2
MPN42103	2.0 mm	2



Solid Core Drill Bits

Item No.	Description	Qty.
MPN11016	1.6 mm	2
MPN11018	1.8 mm	2
MPN11019	2.0 mm	2



AO Style Depth Gauge

Item No.	Description	Qty.
MPN41006	40 mm	1



Ratchet AO/QC Handle

Item No.	Description	Qty.
MSN91001	Cannulated Mini	1



Pickups

Item No.	Description	Qty.
MSN90003	Pickups	1



Retaining Drivers

Item No.	Description	Qty.
MPN30006	T8	4

Section 5



Bone Fragment Pick

Item No.	Description	Qty.
MPN50005	Bone Fragment Pick	1



Retractors

Item No.	Description	Qty.
MPN51006	Mini-Hohmann	2



Compression Slot Drill Guides

Item No.	Description	Qty.
MPN41101	1.6 mm	1
MPN41102	1.8 mm	1
MPN41103	2.0 mm	1



Reduction Forceps

Item No.	Description	Qty.
MPN51015	Mini Pointed	2
MPN50003	Serrated	2



Polyaxial Cone Drill Guide


Item No.	Description	Qty.
MPN42104	1.8/2.0 mm	1

Section 6



Periosteal Elevator

Item No.	Description	Qty.
MPN50006	Periosteal Elevator	1



Bending Pliers

Item No.	Description	Qty.
MPN51008	Bending Pliers	2



5th Met Hook Plate Impactor

Item No.	Description	Qty.
MPN51005	Impactor	1



Plate Cutter

Item No.	Description	Qty.
MPN51007	Plate Cutter	1