


Ancillary Foot Recon Plating System

Tray Layout



UNITE®
FOOT & ANKLE



Section 1

Midfoot/Hindfoot Plates

Item No.	Description	Qty.	Item No.	Description	Qty.
MPP8201L	NC Fusion, Left	1	MPP8501L	Medial Column, Standard, Left	1
MPP8201R	NC Fusion, Right	1	MPP8502L	Medial Column, Long, Left	1
MPP8301U	Middle Column, Standard	1	MPP8501R	Medial Column, Standard, Right	1
MPP8302U	Middle Column, Long	1	MPP8502R	Medial Column, Long, Right	1
MPP8701U	TNC Extension Revision, Universal	1	MPP8601L	Medial Column Extension, Standard, Left	1
MPP8801U	TNCM Extension Revision, Universal	1	MPP8602L	Medial Column Extension, Long, Left	1
MPP5000U	Cotton Osteotomy, Universal	1	MPP8601R	Medial Column Extension, Standard, Right	1
MPP6000U	Evans Osteotomy, Universal	1	MPP8602R	Medial Column Extension, Long, Right	1

Ancillary Foot Recon Plating System

Tray Layout

Section 2

Lisfranc Plates

Item No.	Description	Qty.
MPP9101L	Lisfranc, 1-2 TMT, Tab, Small, Left	1
MPP9102L	Lisfranc, 1-2 TMT, Tab, Medium, Left	1
MPP9103L	Lisfranc, 1-2 TMT, Tab, Large, Left	1
MPP9101R	Lisfranc, 1-2 TMT, Tab, Small, Right	1
MPP9102R	Lisfranc, 1-2 TMT, Tab, Medium, Right	1
MPP9103R	Lisfranc, 1-2 TMT, Tab, Large, Right	1
MPP9001L	Lisfranc, 1-2 TMT, Standard, Small, Left	1
MPP9002L	Lisfranc, 1-2 TMT, Standard, Medium, Left	1
MPP9003L	Lisfranc, 1-2 TMT, Standard, Large, Left	1
MPP9001R	Lisfranc, 1-2 TMT, Standard, Small, Right	1
MPP9002R	Lisfranc, 1-2 TMT, Standard, Medium, Right	1

Item No.	Description	Qty.
MPP9003R	Lisfranc, 1-2 TMT, Standard, Large, Right	1
MPP9201L	Lisfranc, 2-3 TMT, Small, Left	1
MPP9202L	Lisfranc, 2-3 TMT, Medium, Left	1
MPP9203L	Lisfranc, 2-3 TMT, Large, Left	1
MPP9201R	Lisfranc, 2-3 TMT, Small, Right	1
MPP9202R	Lisfranc, 2-3 TMT, Medium, Right	1
MPP9203R	Lisfranc, 2-3 TMT, Large, Right	1
MPP8400L	Lisfranc, 1st TMT with ICJ Arm, Left	1
MPP8400R	Lisfranc, 1st TMT with ICJ Arm, Right	1
MPP8401L	Lisfranc, 2nd TMT with ICJ Arm, Left	1
MPP8401R	Lisfranc, 2nd TMT with ICJ Arm, Right	1

Section 3

Lisfranc Screws

Item No.	Description	Qty.
MSDL3720	Lisfranc Screw, 3.7 X 20mm	2
MSDL3722	Lisfranc Screw, 3.7 X 22mm	2
MSDL3724	Lisfranc Screw, 3.7 X 24mm	2
MSDL3726	Lisfranc Screw, 3.7 X 26mm	2
MSDL3728	Lisfranc Screw, 3.7 X 28mm	2
MSDL3730	Lisfranc Screw, 3.7 X 30mm	2
MSDL3732	Lisfranc Screw, 3.7 X 32mm	2
MSDL3734	Lisfranc Screw, 3.7 X 34mm	2
MSDL3736	Lisfranc Screw, 3.7 X 36mm	2
MSDL3738	Lisfranc Screw, 3.7 X 38mm	2
MSDL3740	Lisfranc Screw, 3.7 X 40mm	2
MSDL3742	Lisfranc Screw, 3.7 X 42mm	2
MSDL3744	Lisfranc Screw, 3.7 X 44mm	2
MSDL3746	Lisfranc Screw, 3.7 X 46mm	2
MSDL3748	Lisfranc Screw, 3.7 X 48mm	2
MSDL3750	Lisfranc Screw, 3.7 X 50mm	2
MSDL3755	Lisfranc Screw, 3.7 X 55mm	2
MSDL3760	Lisfranc Screw, 3.7 X 60mm	2
MSDL4120	Lisfranc Screw, 4.1 X 20mm	2
MSDL4122	Lisfranc Screw, 4.1 X 22mm	2
MSDL4124	Lisfranc Screw, 4.1 X 24mm	2
MSDL4126	Lisfranc Screw, 4.1 X 26mm	2
MSDL4128	Lisfranc Screw, 4.1 X 28mm	2

Item No.	Description	Qty.
MSDL4130	Lisfranc Screw, 4.1 X 30mm	2
MSDL4132	Lisfranc Screw, 4.1 X 32mm	2
MSDL4134	Lisfranc Screw, 4.1 X 34mm	2
MSDL4136	Lisfranc Screw, 4.1 X 36mm	2
MSDL4138	Lisfranc Screw, 4.1 X 38mm	2
MSDL4140	Lisfranc Screw, 4.1 X 40mm	2
MSDL4142	Lisfranc Screw, 4.1 X 42mm	2
MSDL4144	Lisfranc Screw, 4.1 X 44mm	2
MSDL4146	Lisfranc Screw, 4.1 X 46mm	2
MSDL4148	Lisfranc Screw, 4.1 X 48mm	2
MSDL4150	Lisfranc Screw, 4.1 X 50mm	2
MSDL4155	Lisfranc Screw, 4.1 X 55mm	2
MSDL4160	Lisfranc Screw, 4.1 X 60mm	2
MSW04045	Washer	2
MSN10015	Ø2.6mm Cannulated Drill Bit	2
MSN10016	Ø3.7mm Cannulated Over Drill	2
MSN10017	Ø2.8mm Cannulated Drill Bit	2
MSN10018	Ø4.1mm Cannulated Over Drill	2
MSG14150	Ø1.4 x 150mm Guidewire	6
MSN20014	Lisfranc Screw Countersink	1
MPN53002	Lisfranc Tissue Protector	1
MPN53003	Guidewire Sleeve	1
MSN40001	Cannulated Depth Gauge	1

Ancillary Foot Recon Plating System

Tray Layout

Section 4


Plate Screws

Item No.	Description	Qty.
MPSN3512	Screw, Non-Locking, 3.5 X 12mm	2
MPSN3514	Screw, Non-Locking, 3.5 X 14mm	2
MPSN3516	Screw, Non-Locking, 3.5 X 16mm	2
MPSN3518	Screw, Non-Locking, 3.5 X 18mm	2
MPSN3520	Screw, Non-Locking, 3.5 X 20mm	2
MPSN3522	Screw, Non-Locking, 3.5 X 22mm	2
MPSN3524	Screw, Non-Locking, 3.5 X 24mm	2
MPSN3526	Screw, Non-Locking, 3.5 X 26mm	2
MPSN3528	Screw, Non-Locking, 3.5 X 28mm	2
MPSN3530	Screw, Non-Locking, 3.5 X 30mm	2
MPSL3512	Screw, Polyaxial Locking, 3.5 X 12mm	2
MPSL3514	Screw, Polyaxial Locking, 3.5 X 14mm	2
MPSL3516	Screw, Polyaxial Locking, 3.5 X 16mm	2
MPSL3518	Screw, Polyaxial Locking, 3.5 X 18mm	2
MPSL3520	Screw, Polyaxial Locking, 3.5 X 20mm	2
MPSL3522	Screw, Polyaxial Locking, 3.5 X 22mm	2
MPSL3524	Screw, Polyaxial Locking, 3.5 X 24mm	2
MPSL3526	Screw, Polyaxial Locking, 3.5 X 26mm	2
MPSL3528	Screw, Polyaxial Locking, 3.5 X 28mm	2
MPSL3530	Screw, Polyaxial Locking, 3.5 X 30mm	2
MPSN4014	Screw, Non-Locking, 4.0 X 14mm	4
MPSN4016	Screw, Non-Locking, 4.0 X 16mm	4
MPSN4018	Screw, Non-Locking, 4.0 X 18mm	4
MPSN4020	Screw, Non-Locking, 4.0 X 20mm	4
MPSN4022	Screw, Non-Locking, 4.0 X 22mm	4
MPSN4024	Screw, Non-Locking, 4.0 X 24mm	4
MPSN4026	Screw, Non-Locking, 4.0 X 26mm	4
MPSN4028	Screw, Non-Locking, 4.0 X 28mm	4
MPSN4030	Screw, Non-Locking, 4.0 X 30mm	4
MPSN4032	Screw, Non-Locking, 4.0 X 32mm	4
MPSN4034	Screw, Non-Locking, 4.0 X 34mm	4
MPSN4036	Screw, Non-Locking, 4.0 X 36mm	4
MPSN4038	Screw, Non-Locking, 4.0 X 38mm	4
MPSN4040	Screw, Non-Locking, 4.0 X 40mm	4

Item No.	Description	Qty.
MPSN4042	Screw, Non-Locking, 4.0 X 42mm	2
MPSN4044	Screw, Non-Locking, 4.0 X 44mm	2
MPSN4046	Screw, Non-Locking, 4.0 X 46mm	2
MPSN4048	Screw, Non-Locking, 4.0 X 48mm	2
MPSN4050	Screw, Non-Locking, 4.0 X 50mm	2
MPSN4052	Screw, Non-Locking, 4.0 X 52mm	2
MPSN4054	Screw, Non-Locking, 4.0 X 54mm	2
MPSN4056	Screw, Non-Locking, 4.0 X 56mm	2
MPSN4058	Screw, Non-Locking, 4.0 X 58mm	2
MPSN4060	Screw, Non-Locking, 4.0 X 60mm	2
MPSL4014	Screw, Polyaxial Locking, 4.0 X 14mm	4
MPSL4016	Screw, Polyaxial Locking, 4.0 X 16mm	4
MPSL4018	Screw, Polyaxial Locking, 4.0 X 18mm	4
MPSL4020	Screw, Polyaxial Locking, 4.0 X 20mm	4
MPSL4022	Screw, Polyaxial Locking, 4.0 X 22mm	4
MPSL4024	Screw, Polyaxial Locking, 4.0 X 24mm	4
MPSL4026	Screw, Polyaxial Locking, 4.0 X 26mm	4
MPSL4028	Screw, Polyaxial Locking, 4.0 X 28mm	4
MPSL4030	Screw, Polyaxial Locking, 4.0 X 30mm	4
MPSL4032	Screw, Polyaxial Locking, 4.0 X 32mm	4
MPSL4034	Screw, Polyaxial Locking, 4.0 X 34mm	4
MPSL4036	Screw, Polyaxial Locking, 4.0 X 36mm	4
MPSL4038	Screw, Polyaxial Locking, 4.0 X 38mm	4
MPSL4040	Screw, Polyaxial Locking, 4.0 X 40mm	4
MPSL4042	Screw, Polyaxial Locking, 4.0 X 42mm	2
MPSL4044	Screw, Polyaxial Locking, 4.0 X 44mm	2
MPSL4046	Screw, Polyaxial Locking, 4.0 X 46mm	2
MPSL4048	Screw, Polyaxial Locking, 4.0 X 48mm	2
MPSL4050	Screw, Polyaxial Locking, 4.0 X 50mm	2
MPSL4055	Screw, Polyaxial Locking, 4.0 X 55mm	2
MPSL4060	Screw, Polyaxial Locking, 4.0 X 60mm	2
MPN10030	Ø3.0mm Drill Bit	2
MPN40008	Ø3.0mm Locking Drill Guide Tower	2
MPN40009	Ø3.0mm Compression Slot Drill Guide	1

Ancillary Foot Recon Plating System

Tray Layout



Section 5

Instrumentation Tray

Item No.	Description	Qty.
MPN51006	Mini-Hohmann Retractor	2
MPN50005	Bone Fragment Pick	1
MPN50006	Periosteal Elevator	1
MPN53000	Pointed Reduction Clamp, Long Throw	1
MPN50022	Pointed Reduction Clamp, Wheel	1
MPN53001	Lisfranc Targeting Guide Clamp	1

To schedule a case, contact your Medline UNITE Representative or visit medlineunite.com for more information.