


Distal Tibia Plating System

Tray Layout—Level 1



UNITE
FOOT & ANKLE



Section 1

Anterior Distal Tibia Plates

Item No.	Description	Qty.
MPDT001U	3 Hole	1
MPDT002U	6 Hole	1
MPDT003U	9 Hole	1

Posterolateral Distal Tibia Plates

Item No.	Description	Qty.
MPDT301L	3 Hole, Left	1
MPDT301R	3 Hole, Right	1
MPDT302L	5 Hole, Left	1
MPDT302R	5 Hole, Right	1
MPDT303L	7 Hole, Left	1
MPDT303R	7 Hole, Right	1
MPDT304L	9 Hole, Left	1
MPDT304R	9 Hole, Right	1

Section 2

Medial Distal Tibia Plates

Item No.	Description	Qty.
MPDT201L	4 Hole, Left	1
MPDT201R	4 Hole, Right	1
MPDT202L	6 Hole, Left	1
MPDT202R	6 Hole, Right	1
MPDT203L	8 Hole, Left	1
MPDT203R	8 Hole, Right	1
MPDT204L	10 Hole, Left	1
MPDT204R	10 Hole, Right	1
MPDT205L	12 Hole, Left	1
MPDT205R	12 Hole, Right	1
MPDT206L	14 Hole, Left	1
MPDT206R	14 Hole, Right	1


Section 3

Anterolateral Distal Tibia Plates

Item No.	Description	Qty.
MPDT101L	4 Hole, Left	1
MPDT102L	6 Hole, Left	1
MPDT103L	8 Hole, Left	1
MPDT104L	10 Hole, Left	1
MPDT105L	12 Hole, Left	1
MPDT106L	14 Hole, Left	1
MPDT101R	4 Hole, Right	1
MPDT102R	6 Hole, Right	1
MPDT103R	8 Hole, Right	1
MPDT104R	10 Hole, Right	1
MPDT105R	12 Hole, Right	1
MPDT106R	14 Hole, Right	1

Distal Tibia Plating System

Tray Layout—Level 1



Section 4

3.5 mm Polyaxial Locking Screws

Item No.	Description	Qty.
MPSL3522	3.5 x 22 mm	2
MPSL3524	3.5 x 24 mm	2
MPSL3526	3.5 x 26 mm	2
MPSL3528	3.5 x 28 mm	2
MPSL3530	3.5 x 30 mm	4
MPSL3532	3.5 x 32 mm	4
MPSL3534	3.5 x 34 mm	4
MPSL3536	3.5 x 36 mm	4
MPSL3538	3.5 x 38 mm	4
MPSL3540	3.5 x 40 mm	4
MPSL3542	3.5 x 42 mm	4
MPSL3544	3.5 x 44 mm	4
MPSL3546	3.5 x 46 mm	4
MPSL3548	3.5 x 48 mm	4
MPSL3550	3.5 x 50 mm	4
MPSL3555	3.5 x 55 mm	2
MPSL3560	3.5 x 60 mm	2

3.5 mm Non-Locking Screws

Item No.	Description	Qty.
MPSN3522	3.5 x 22 mm	2
MPSN3524	3.5 x 24 mm	2
MPSN3526	3.5 x 26 mm	2
MPSN3528	3.5 x 28 mm	2
MPSN3530	3.5 x 30 mm	4
MPSN3532	3.5 x 32 mm	4
MPSN3534	3.5 x 34 mm	4
MPSN3536	3.5 x 36 mm	4
MPSN3538	3.5 x 38 mm	4
MPSN3540	3.5 x 40 mm	4
MPSN3542	3.5 x 42 mm	4
MPSN3544	3.5 x 44 mm	4
MPSN3546	3.5 x 46 mm	4
MPSN3548	3.5 x 48 mm	4
MPSN3550	3.5 x 50 mm	4
MPSN3555	3.5 x 55 mm	2
MPSN3560	3.5 x 60 mm	2

3.7/4.0 mm IM Fibula Implant


Item No.	Description	Qty.
MSFB0065	3.7/4.0 x 65 mm	1
MSFB0070	3.7/4.0 x 70 mm	1
MSFB0075	3.7/4.0 x 75 mm	1
MSFB0080	3.7/4.0 x 80 mm	1
MSFB0085	3.7/4.0 x 85 mm	1
MSFB0090	3.7/4.0 x 90 mm	1
MSFB0095	3.7/4.0 x 95 mm	1
MSFB0100	3.7/4.0 x 100 mm	1
MSFB0110	3.7/4.0 x 110 mm	1
MSFB0120	3.7/4.0 x 120 mm	1
MSFB0130	3.7/4.0 x 130 mm	1
MSFB0140	3.7/4.0 x 140 mm	1
MSFB0150	3.7/4.0 x 150 mm	1

Instruments

Item No.	Description	Qty.
MPN50019	Ball Spike Pusher, AO/QC	1

Distal Tibia Plating System

Tray Layout—Level 2



Section 5

Item No.	Description	Qty.
MPN50022	Forceps, Reduction Forceps, Pointed, Wheel	1
MSN90006	Non-Ratcheting Cannulated Handle, Standard	1
MPN50021	Distal Tibia Plates Inserter	1
MPN50029	Curved Banana Elevator, AO/QC	1
MPN40002	2.8 mm Locking Tower Drill Guide	2

Distal Tibia XL Plates Caddy*

Item No.	Description	Qty.
MPDT107L	16 Hole Anterolateral Distal Tibia Plate, Left	1
MPDT107R	16 Hole Anterolateral Distal Tibia Plate, Right	1
MPDT108L	18 Hole Anterolateral Distal Tibia Plate, Left	1
MPDT108R	18 Hole Anterolateral Distal Tibia Plate, Right	1
MPDT207L	16 Hole Medial Distal Tibia Plate, Left	1
MPDT207R	16 Hole Medial Distal Tibia Plate, Right	1
MPDT208L	18 Hole Medial Distal Tibia Plate, Left	1
MPDT208R	18 Hole Medial Distal Tibia Plate, Right	1

**This is an ancillary caddy that is not standard in the Distal Tibia Plating System and is available upon request. Caddy is not pictured.*

To schedule a case, contact your Medline UNITE Representative or visit medlineunite.com for more information.