

ANCHORED IN PRECISION


FORCEWEB™ Synthetic Ligament Augmentation

Constructed from an open weave synthetic, non-resorbable material, the implant's dense microfilament structure allows for wicking of biological fluids to facilitate cellular adherence and tissue in-growth.¹






FORCECORD™ Midsubstance Repair System

This densely woven tubular implant outperforms 2 mm suture tape in both ultimate tensile strength and pull-out strength.^{1,2,3}



DEXLOCK® Insertional Achilles Repair Kit

Sterile convenience kit includes four Ø4.5 mm DEXLOCK Knotless Anchors, two of which are preloaded with 2 mm tape and needles, and all necessary instrumentation.



UNITE®
FOOT & ANKLE



Explore our product portfolio by
system, anatomy or procedure.

REFERENCE: 1. Qual-148371, 2. Qual-131707, 3. Qual-143158 / © 2026 Medline. All rights reserved. Medline, Medline UNITE, and DEXLOCK are registered trademarks, FORCEWEB and FORCECORD are trademarks of Medline Industries, LP. / MKT19W9038758 / e25602 / 49